



LOUPES

Magnification Systems

SIZE

66□15
Western Fitting Japan/Asian Fitting

ACCESSORIES INCLUDED

Elastic band

FEATURES

Quick Loupes are loupes accessible to everyone. Easy to use with predefined interpupillary distances that can be selected according to your needs, they maintain the quality, comfort and optical characteristics of a Tailor Made loupes.
Available only in Galilean 2.5X version.

Quick Loupes ist die Lupenbrille, die für alle erschwinglich ist. Einfach zu bedienen, mit vordefinierten Pupillenabständen, die Sie nach Ihren Bedürfnissen wählen können, gewährleistet es alle Vorteile, den Komfort und die optischen Eigenschaften der maßgeschneiderten Geräten.
Nur in der galileische 2.5X Version verfügbar

INTERPUPILLARY DISTANCE

Size	PD at infinite viewing (mm)	Range covered (mm)
S	60	58 - 62
M	64	62 - 66
L	68	66 - 70

COLOURS CHART

available colour combinations for frame/optics

Black Green




White Pink



Quick Loupes

White Red



Quick
Available Galilean 2.5X

www.univetloupes.com



UNIVET LOUPES Spa Via Giovanni Prati, 87 - 25086 Rezzato (BS) Italy +39 030 2499411 · info@univetloupes.com